

Thermal Therapy acts to re-establish cell metabolism balance, making tissues more compact and younger looking. It improves microcirculation, reduces oedema, fat accumulation and skin laxity.

Thermal Therapy provides a painless non invasive treatment with surprising results from the first session.

TECHNICAL CHARACTERISTICS

- working frequencies: 420-500-720-1000 khz to optimize the treatment to various tissues and at various depths
- Measurement and transferred energy indication to the tissues (J) for a correct and repeatable energy dosage
- Continuous and pulsed modes of operation (for athermal treatments)
- Power from 80 up to 300 w, continuously adjustable for different options of use
- High-quality touch screen colour display and power control knob
- Up to 4 handpieces: bipolar face handpiece, bipolar body handpiece, monopolar face handpiece and monopolar body handpiece



REG. N. 1889
UNI EN ISO 9001: 2008
UNI EN ISO 13485: 2012



GP
General Project
ADVANCED BIOMEDICAL TECHNOLOGIES

General Project s.r.l.

Via della Gora 15/19
50025 Montespertoli (Firenze) Italy
Tel: +39 0571 675076
Fax: +39 0571 675077
www.generalproject.com
info@generalproject.com

82022630UK rev.0



Revolutionary Heat

Thermal Therapy



thermal therapy

Customizable face & body Thermal Treatments with a wide range of choices for a large number of treatment options

Up to 4 handpieces, for face and body, allow all thermal non invasive applications which exploit the diathermy principle, the deep heating of the skin tissue, due to high frequency electromagnetic waves. It is meant to treat skin laxity, wrinkles and face and body flaws. Its lifting effect is produced by the temperature increase which causes a contractions of the collagen fibres with a reshaping of the treated area. Besides, the vascularisation and liquid drainage effects can be focused on the most superficial layers, thus improving skin tone and glow and reducing localized fat and peau d'orange effect to obtain a brighter and smoother skin with visible results from the first session.

The Smart model works with 2 bipolar RF handpieces, one for the face and the other for the body. The generated heat causes the local contraction of collagen fibres and improves the skin's laxity and appearance, remodelling and enhancing dermal density.

The unit can be completed with the monopolar handpieces (capacitive and resistive), which, thanks to the depth of their action, allow the heating of the dermis deepest tissues and of the fat tissues. They are particularly advised to treat cellulite, to perform body shaping and deep tissue regeneration.

Furthermore, the possibility to select among **4 treatment frequencies** for each handpiece allow to vary the action depth depending on the blemish or the area to treat.

In pulse mode, the unit works also in diathermy.

A wide range of combinations and relevant treatments make Thermal Therapy a device able to fulfill all possible needs.

APPLICATIONS

- Wrinkles reduction
- Firming, Lifting and rejuvenation effect
- Facial contour redefinition
- Body shaping and lifting



Face Handpiece

Body Handpiece

